



Innovative
Healthcare made Accessible

—
PRODUCT CATALOGUE



OUR MOTIVATION

At Drive DeVilbiss Healthcare we are motivated to work together with our customers and partners, in order to make everyday life easier, for our product end users. That is why we give our very best, every day to insure we are supplying innovative products and excellent customer service. Our friendly and knowledgeable employees look forward to supporting you.

OUR STORY

Drive was formed in January 2000 under the name Medical Depot, Inc. In 2015 DeVilbiss Healthcare was acquired by Drive and the company became known as Drive DeVilbiss Healthcare. The continuous growth and expansion throughout the years is supported by the best professionals in the business, state-of-the-art technology and a tenacious winning attitude. This culture and drive lead to the significant growth and formation of a multinational organization.

EXCELLENT QUALITY
WORLD CLASS PRODUCT

OVER 150 HIGHLY MOTIVATED,
KNOWLEDGEABLE EMPLOYEES

MORE THAN 24
LOCATIONS WORLDWIDE

DRIVE DEVILBISS INTERNATIONAL
SUPPLY MORE THAN 60 COUNTRIES
WORLDWIDE



MOBILITY



ROLLATOR



BATHROOM



SCOOTER



CANES



Multifunctional wheelchair Multitec

MOBILITY

- Adjustable-angle seat unit continuously supported by gas pressure damper
- Adjustable-angle, ergonomically shaped backrest continuously supported by gas pressure damper independently of the seat unit
- Specially shaped surface structure, optimal pressure distribution, and stable positioning of the user
- Height-adjustable push handles
- Folding leg supports with knee and calf cushions (height and depth adjustable)
- Anti-tilt supports (height-adjustable) and headrest (optimally adjustable height and inclination) fitted as standard
- Puncture-proof 24" PU drive wheels with quick-release axles and chrome-plated hand rims
- Accessories: abductor wedge, side supports, wheelchair tray, adapter plates to widen seat, and adapter set for pushing aid PowerStroll
- Care aid number: 50,45,07,2033
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems



Adjustable headrest



Height-adjustable push handles



Height-adjustable tilt aid



Product	Multifunctional wheelchair Multitec		
Seat width	39 cm	44 cm	49 cm
without drum brake	900950100	900950200	900950300
with drum brake	900900100	900900200	900900300
Overall width without / with DB	60 cm / 65 cm	65 cm / 70 cm	70 cm / 80 cm
Weight without / with DB	31 kg / 32,5 kg	32 kg / 33,5 kg	34,5 kg / 13 kg

Technical Data	
Overall length	112 - 152 cm
Overall height	97 - 114 cm
Seat depth	42 - 48 cm
Seat height	45 cm
Seat angle from	-1° bis 19°
Backrest height	51 - 65 cm
Backrest angle	-3° bis 36°
Footrest adjustment	40 - 60 cm
Max. Load	135 kg

Standard wheelchair Ecotec 2G

MOBILITY 

- Very light and versatile wheelchair
- Armrests can be swiveled away
- Folding footrests with heel straps
- Upholstered seat back and seat cushion
- Height of aluminum wheel forks can be adjusted to three positions
- Powder-coated steel hand rims
- Seat height adjustment of 46 and 51 cm
- Quick-release axles as standard for fast wheel installation or removal
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Integrated brake lever extension can be extended and folded over
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems

WITH REFINED IDEAS



CRASHTEST
ISO 7176-19

Wheel fork



Adjustable back-rest



Product	Standard wheelchair Ecotec 2G			
Seat width	38 cm	42 cm	46 cm	50 cm
without drum brake	910300100	910300200	910300300	910300400
with drum brake	910350100	910350200	910350300	910350400
Overall width without / with DB	56,5 cm / 59 cm	60,5 cm / 63 cm	64,5 cm / 67 cm	68,5 cm / 71 cm
Width folded without / with DB	30 cm / 32 cm	30 cm / 32 cm	30 cm / 32 cm	30 cm / 32 cm
Weight without / with DB	18,6 kg / 19,6 kg	18,8 kg / 19,8 kg	19 kg / 20 kg	19,2 kg / 20,2 kg

Technical Data Allgemein	
Overall length o. / m. TB	92 cm / 107 cm
Overall height	90 - 95 cm
Seat depth	42 cm
Seat height	49 / 54 cm
Backrest height	39 cm
Footrest adjustment	36 - 47 cm
Max. Load	130 kg

Lightweight wheelchair Litec 2G

MOBILITY 

- Very light and versatile wheelchair
- Armrests can be swiveled away
- Adjustable-angle (0° - 30°) and folding footrests as standard
- Leg support height adjustable to six positions
- Upholstered seat back and seat cushion
- Reinforced aluminum frame for optimal stability
- Height of aluminum wheel forks can be adjusted to three positions
- Three-position seat height adjustment from 46 to 52 cm
- Wheelbase extension possible
- 175 x 45 mm PU steering wheels and 24" full rubber drive wheels
- Quick-release axles as standard for fast wheel installation or removal
- Crash test according to ISO 7176-19: approved for passenger transport in vehicles in combination with all approved passenger and wheelchair restraint systems

A WHOLE RANGE OF OPTIONS FOR
INDIVIDUAL ADJUSTMENT



Armrest can be swiveled



Wheel fork



Foot plate with heel strap



Product	Standard wheelchair Ecotec 2G			
Seat width	38 cm	42 cm	46 cm	50 cm
without drum brake	910300100	910300200	910300300	910300400
with drum brake	910350100	910350200	910350300	910350400
Overall width without / with DB	56,5 cm / 59 cm	60,5 cm / 63 cm	64,5 cm / 67 cm	68,5 cm / 71 cm
Width folded without / with DB	30 cm / 32 cm	30 cm / 32 cm	30 cm / 32 cm	30 cm / 32 cm
Weight without / with DB	18,6 kg / 19,6 kg	18,8 kg / 19,8 kg	19 kg / 20 kg	19,2 kg / 20,2 kg

Technical Data Allgemein	
Overall length o. / m. TB	92 cm / 107 cm
Overall height	90 - 95 cm
Seat depth	42 cm
Seat height	49 / 54 cm
Backrest height	39 cm
Footrest adjustment	36 - 47 cm
Max. Load	130 kg

Lightweight Rollator Gigo 2G

ROLLATOR

- Particularly light and attractive aluminum rollator
- One-hand folding mechanism ensures safe and comfortable folding
- Ergonomic handles can be adjusted to 6 height positions
- Max. load of basket and tray 5 kg each
- Supplied as standard with basket, tray, backrest, and cane holder



Double-function brake

Backrest



Product	Gigo 2G
Order number	735100200

Technical Data	
Overall length x Overall width	61 cm x 61,5 cm
Handle height	87 - 96 cm
Seat width x Seat depth	34 cm x 33,5 cm
Seat height	54 cm
Width between handles	40,5 cm
Folding Dimensions	28 x 61,5 x 100 cm
Weight	6,6 kg
Max. Load	120 kg

Rollator Cristallo²

ROLLATOR

- For rheumatism and arthritis patients
- Cushioned and individually adjustable forearm supports
- Automatic lock to prevent folding
- Multiple reflectors ensure safety
- Supplied as standard with basket, tray, and cane holder



Forearm support, padded

Forearm support with adjustable length



Product	Cristallo ²
Order number	722500000

Technical Data	
Overall length x Overall width	67 cm x 65 cm
Armrest height	102 - 118 cm
Seat width x Seat depth	38 cm x 16 cm
Seat height	61 cm
Width between handles	38 cm
Folding Dimensions	23 x 58 x 74 cm
Weight	12,3 kg
Max. Load	120 kg

Indoor-Rollator Roomba

ROLLATOR 

The sporty companion for use at home

- Designed specifically for use indoors
- Can be used anywhere thanks to its narrow design, even when space is limited
- Soft hand grips with end-to-end brake lever
- Brake can be operated and locked with one hand
- Supplied as standard with basket and large tray with round depression



Large tray



Versatile mesh bag for personal items



Safety latch



Product		Roomba
Order number		730200100
Technical Data		
Overall length x Overall width		51,5 cm x 57 cm
Handle height		81,4 - 94,1 cm
Tray width		45 cm
Tray depth		30 cm
Width between rear wheels		44 cm
Folding Dimensions		29 x 55 x 83 cm
Weight		5,6 kg
Max. Load		100 kg

Rollator Nitro

ROLLATOR

- Can be used indoors and outdoors
- Increased stability and safety thanks to double cross struts
- Extra-large 25 cm steering wheels for optimal steering and pushing comfort
- Passive lighting at the rear, front, and sides as standard
- Cushioned seat surface and backrest as standard
- Simple handle height adjustment
- Cane holder included as standard
- Ergonomic, height-adjustable handles
- Dual-function brake (for braking and parking)



SIZE M ALSO IN RED



THE NITRO FAMILY AVAILABLE IN 4 COLOURS AND 4 SIZES



Product				Size L
champagner	724100300	724100000	724200000	724300000
white	-	-	724200100	724300100
black	-	-	724200200	-
red	-	-	724200400	-

Technical Data				
Overall length x Overall width	70,5 cm x 52 cm	70,5 cm x 61 cm	70,5 cm x 61 cm	70,5 cm x 61 cm
Handle height	71,8 - 77 cm	79 - 91 cm	84 - 96,5 cm	89 - 101 cm
Seat width x Seat depth	37 cm x 20 cm	44,5 cm x 22 cm	44,5 cm x 22 cm	44,5 cm x 22 cm
Seat height	45 cm	52 cm	57 cm	62 cm
Width between handles	34,5 cm	45 cm	45 cm	45 cm
Folding Dimensions	70,5 x 25 x 71,8 cm	70,5 x 25 x 82 cm	70,5 x 25 x 86 cm	70,5 x 25 x 92 cm
Total weight	6,7 kg	7,1 kg	7,2 kg	7,3 kg
Max. Load	130 kg	130 kg	130 kg	130 kg

Shower Stool Delphi Plus

BATHROOM

- The hygiene opening makes daily body care easier
- Safe, stable, and space-saving
- Seat surface with drainage slits
- Low weight
- Three-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Anodized aluminum tube legs
- Easy to clean
- Delivery in practical take-away packaging

Product		Delphi Plus
Order number		540200000
Technical Data		
Overall width		48 cm
Overall length		48 cm
Overall height		55 cm
Seat height		45 - 55 cm
Seat plate diameter		38 cm
Total weight		1,6 kg
Max. Load		180 kg



Shower Stool DS 130

BATHROOM

- Comfortable, safe, and stable
- Large seat surface with drainage holes
- Hygienic cutaway simplifies personal hygiene
- Five-position seat height adjustment through telescopic legs
- Stable thanks to rubber feet
- Powder-coated aluminum frame
- Easy to clean

Hygiene opening Height-adjustable legs



Product		DS 130
Order number		550200200
Technical Data		
Overall width min. / max.		53 cm / 54 cm
Overall height min. / max.		63 cm / 72 cm
Overall depth min. / max.		55 cm / 57 cm
Seat width		41 cm
Seat depth		41 cm
Seat height		42 - 52 cm
Total weight		2,8 kg
Max. Load		130 kg



Shower Chair DSF 130

BATHROOM

- Easy to stow - folding
- Space-saving thanks to folding mechanism
- Large seat surface with drainage holes
- Hygienic cutaway simplifies personal hygiene
- Comfortable, safe, and stable
- Five-position seat height adjustment through telescopic legs



Product	DSF 130
Order number	550200100

Technical Data	
Overall width minimal / maximal	52 cm / 53 cm
Overall height minimal / maximal	90 cm / 100 cm
Overall depth minimal / maximal	56 cm / 58 cm
Seat width x Seat depth	41 cm x 41 cm
Seat height	42 - 52 cm
Backrest height	48 cm
Total weight	3,8 kg
Max. Load	130 kg



Shower Chair Duro KD

BATHROOM

- Comfortable, safe, stable, can be dismantled
- Large seat surface with drainage holes and carrying handles
- Low weight
- Available with and without backrest
- Five-position seat height adjustment through telescopic legs
- Anodized aluminum tube legs



Product	with backrest	
Order number	550500900	550500950

Technical Data		
Overall width	53 cm	53 cm
Overall height	54 cm	81 cm
Overall depth	52 cm	56 cm
Seat width	50 cm	50 cm
Seat depth	30 cm	30 cm
Seat height	42 - 52 cm	42 - 52 cm
Total weight	2 kg	2,9 kg
Max. Load	130 kg	130 kg



Raised Toilet Seat TSE 120 Plus

BATHROOM

- NEW: with innovative top mounting
- Two hygienic cutaways and integrated seat
- Recessed grip on the lid makes it easy to open
- Folding, non-slip armrests
- Simple, tool-free assembly
- Three-position height adjustment
- Smooth surfaces simplify cleaning
- Delivery includes assembly gloves
- Delivery in practical take-away packaging



Product	TSE 120 Plus
Order number	520500500

Technical Data	
Overall width	55 cm
Overall length	49 cm
Height adjustment	8,5 / 12 / 15,5 cm
Width between armrest	46 cm
Opening length	27 cm
Opening width	20,5 cm
Total weight	4,1 kg
Max. Load	120 kg

Raised Toilet Seat Ticco 2G/10

BATHROOM

- Raises the toilet seat by 10 cm
- Ergonomic shape
- Maximum opening with two hygienic cutaways
- Recessed grip on the lid makes it easy to open
- Simple, tool-free assembly
- Fixation clamps enable perfect hold
- Smooth surfaces for optimal cleaning and disinfection
- Lightweight with maximum load capacity
- Delivery in practical take-away packaging



Product	without lid	
Order number	523400000	523300000

Technical Data		
Overall width	39 cm	39 cm
Overall length	39 cm	39 cm
Height front	10 cm	8,6 cm
Height rear	14,6 cm	12 cm
Opening length	28,5 cm	28,5 cm
Opening width	21,5 cm	21,5 cm
Total weight	1,35 kg	1,05 kg
Max. Load	225 kg	225 kg

Alternating pressure system MED AIRE DuoWave

BATHROOM

- Individual adjustment of the cell pressure to the weight of the patient, both lying and sitting
- Individually usable thanks to possible static pressure-relief positioning
- Continuous, manual cell pressure adjustment based on a weight scale
- Bi-elastic, quilted, and breathable cover
- Fast repair of the overlay is possible thanks to individual air cell exchange
- Recommend as therapeutic support for all decubitus categories (1 to 4) and for prophylaxis in patients up to highdecubitus risk



Product		MED AIRE DuoWave
Order number		82000200
Technical Data		
Pump		
Dimensions (WxDxH)		28 x 10 x 20 cm
Total weight		2,2 kg
Power Supply		2v30 V, 50 Hz
Input		12 W
Pressure range		30 - 60 mbar
Mattress		
Dimensions (LxWxH)		200 x 90 x 14 cm
User weight		36 kg - 150 kg



integrated foam pocket in the base

Comfort cover

CPR function



Suction Cup Grab Bar Solido / Solido T

BATHROOM

- The mobile all-rounder for use at home or on the go
- Makes it easier to hold on and stand up in the bathtub or shower
- Versatile and quickly ready to use on all smooth surfaces
- Vacuum rocker lever generates extra-strong suction with a tensile force of 70 kg
- Available in three lengths: 32.5, 58, and 88 cm
- Delivered as complete unit ready for use
- Lieferung in praktischer Mitnahmeverpackung

Vacuum rocker lever

Solido, short

easily adjustable length



Product	SOLIDO			SOLIDO T	
	Length 32.5	Length 58 cm	Length 88 cm	short	medium
Order number	660100100	660100200	660100300	660100400	660100500
Technical Data					
	SOLIDO			SOLIDO T	
Overall length	32,5 cm	58 cm	88 cm	44 - 57 cm	56 - 69 cm
Gripping surface	10 cm	32 cm	62 cm	34 cm	46 cm
Wall clearance	3,5 cm	3,5 cm	3,5 cm	3,5 cm	3,5 cm
Suction Cup Ø	12 cm	12 cm	12 cm	12 cm	12 cm
Total weight	600 g	1 kg	1,3 kg	800 g	900 g
Max. Load	70 kg	70 kg	70 kg	70 kg	70 kg



Raised toilet seat Ticco 2G/15

BATHROOM

- Raises the toilet seat by 15 cm
- Ergonomic shape
- Maximum opening with two hygienic cutaways
- Recessed grip on the lid makes it easy to open
- Simple, tool-free assembly
- Fixation clamps enable perfect hold
- Smooth surfaces for optimal cleaning and disinfection
- Lightweight with maximum load capacity
- Delivery in practical take-away packaging

Product	with lid	without lid
Order number	523600400	523600300

Technical Data		
Overall width	39 cm	39 cm
Overall length	39 cm	39 cm
Height front	15,3 cm	13,6 cm
Height rear	16,6 cm	15,1 cm
Opening length	27,3 cm	27,3 cm
Opening width	20,3 cm	20,3 cm
Total weight	1,71 kg	1,44 kg
Max. Load	225 kg	225 kg



Two hygiene openings



Fixing screw for easy assembly



also available with a lid



Wheeled Commode TRS 200 XXL

BATHROOM

- Wide and stable wheeled commode for maximum requirements
- Parking brake on the rear casters
- Fold-down armrests simplify transfer
- Stable, folding, and removable footrests
- With receptacle as standard

With integrated cover plate



Adjustable footrest



Bucket with odor trap



Product	TRS 200 XXL
Order number	530200000

Technical Data	
Overall width	72 cm
Overall length	96 cm
Overall height	98 cm
Seat width x Seat depth	56 cm x 42 cm
Seat height	52 cm
Footrest adjustment	42 - 57 cm
Clearance height	44 cm
Total weight	15,3 kg
Max. Load	200 kg



Scooter RAVEN (BL270 SCOUT AIR)

SCOOTER

- Can be disassembled quickly without the need for tools
- Low net weight
- Clear control panel
- Small turning circle
- Standard pneumatic tyres
- Anti-tip wheels



Product		BL270 Scout Air
red		MS025RD-G-AIR
Technical Data		
Speed		6 km/h
Battery capacity 2x		20 Ah
Motor capacity / Peak power		270 W / 500 W
Turning radius		99 cm
Gradeability		14 %
Clearance		6 cm
Front tires / hinten		200 x 50 mm
Dimensions (LxBxH)		108 x 48 x 91 cm
Seat width x Seat depth		43 cm x 36 cm
Width between armrest		50 cm
Weight without battery		34 kg
Max. Load		130 kg



Einfaches Bedienteil



Luftbereifung



Einfach zerlegbar



Quad Cane M

CANES

- The quad cane offers maximum safety and ensures optimal stability
- The ergonomic hand grip sits comfortably and securely in the hand
- The height can be individually adjusted to the user
- The grip thrills with its high-quality finish

Secure and optimal footing



Product		Quad cane M
Order number		WS016
Technical Data		
Handle height		66 - 92 cm
Total weight		750 g
Max. Load		115 kg





SCAN FOR MORE DETAILS

Get In Touch With Us

WWW.SANRAIMED.COM

ADDRESS

Office Tower B, Elante, Office No. 506
Plot No.178-178A Industrial Area Phase-1
Chandigarh- 160002

CUSTOMER SUPPORT

P: +91 91 2957 5666
E: customersupport@sanraimed.com