

OxyFlow 5

OXYGEN CONCENTRATOR

- 1-5L/min
- Lightweight and Compact Design
- Economically Priced
- Fast & Easy Maintenance
- Low Purity Alarm



OxyFlow 5

Specification	Model: Oxyflow
Oxygen Flow	0.5 to 5 Liters per minute
Delivery Method	Continuous
Oxygen Purity	93% ± 3% @ 5LPM 93% ± 3% @ 3LPM
Weight (kg)	14.2kg
Dimensions	600H, 390W, 230D(mm)
Outlet Pressure	7.2-7.5 Psi
Power	300 Watts
Noise Level	≤45dB
O2 Purity Alarm	Low Oxygen
Electrical Requirements	110V, 60HZ / 220V, 50HZ

