



FINDING THE RIGHT SLEEP APNEA MASK IS ESSENTIAL FOR EFFECTIVE SLEEP THERAPY.



Nasal mask AirFit N20

- Comfortably moves with you throughout the night, regardless of sleep position
- Features InfinitySeal™ cushion technology, our most adaptive yet
- Ideal for those who breathe through their nose during sleep



Resmed AirFit P30i Nasal Pillow Mask

The ResMed AirFit P30i pillows mask features a top-of-the-head tube design and pillows cushion that makes it easier to sleep in many different positions! Enjoy the freedom to move in your sleep and get closer to your bed partner.

Recommended sizing: small starter pack for females and standard starter pack for males.



Resmed AirFit F30


Small, quiet and designed for freedom, the AirFit F30 CPAP mask delivers full-face functionality in an ultra-compact format. With nothing on your nose, there's no risk of red marks or discomfort on the bridge of your nose. Patients experience lower obstruction on face and can comfortably wear their spectacles with AirFit F30 full-face mask. It is compatible with all ResMed CPAP devices along with AirMini Autoset ravel device.



Resmed AirFit N30i Nasal Cradle Mask

The ResMed AirFit N30i nasal cradle mask sits under your nose and features a top-of-the-head tube design that makes it easier to sleep in any position! Red marks and discomfort are a thing of the past. Enjoy the freedom to move in your sleep and get close to your bed.

AirFit N30i comes in 2 sizes – Standard & Small (each pack contains 3 cushions in it).



Resmed AirFit P10 Nasal Pillow Mask

The ResMed AirFit™ P10 nasal pillows mask is our quietest mask yet. Bringing you the latest technology in sleep apnea therapy, the AirFit P10 is approximately 50% quieter and 50% lighter than its predecessor.

A recent ResMed study found that patients using the AirFit P10 wore it for 40 minutes longer per night. Offering premium comfort and ease of use without compromising stability, the AirFit P10 also sports a unique split-strap headgear design for minimal facial contact and a real sense of freedom.



AcuCare F1-4 hospital vented full face mask

- Vented mask with a built-in AAV for leak ventilation (single limb).
- Colour-coded, easily identifiable, and designed to fit with common breathing circuits used in hospitals around the world.
- Clear elbow indicates mask contains an anti-asphyxia valve (AAV); can be used without a vented circuit.
- Uses 'Full face' mask setting in ResMed ventilators.
- Single patient use for up to 7 days.
- Sizing template located on packaging for quick mask selection.



AirFit™ P10 for AirMini™ Nasal pillows mask

- Compact, under-the-nose frame for a clear field of vision
- Simple design for easy fitting, cleaning and assembly
- HumidX™ waterless humidification designed to provide comfort and convenience
- ActiveAir™ venting for compact, portable therapy
- Designed exclusively for AirMini



Full face mask AirFit F20

- Specially-designed to fit all patients, regardless of facial profile
- Designed to perform across a wide range of therapy pressures
- Modular frame fits all cushion sizes